MACSQuant[®] Instrument

User manual

Version 5
Original instructions

Copyright © 2013 Miltenyi Biotec. All rights reserved.

No part of this user manual may be reproduced, stored in a retrieval system, transmitted, published, or distributed in any form or by any means, electronically, mechanically, by photocopying, microfilming, recording, or otherwise, without the prior written consent of Miltenyi Biotec; however, notwithstanding the foregoing, the owners of the MACSQuant Instrument may make copies solely for purposes of training personnel in the use and servicing of the unit within their business or organization. Maximal care has been taken by Miltenyi Biotec in the preparation of this user manual. However, Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this document. The information in this document is provided "as is" without warranty of any kind, and is subject to change without notice. MACS, MACSQuant, and MACSQuantify are registered trademarks or trademarks of Miltenyi Biotec GmbH. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only. Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use.

MACSQuant[®] Instrument

User manual

Version 5 EN

Original instructions

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de
www.miltenyibiotec.com

Read the user manual before using the instrument

Before using this instrument read the chapter **Important safety information** and all other information contained in this user manual, including all safety and operating instructions. Pay special attention to warnings displayed on the instrument. Failure to read and follow these guidelines could lead to improper or incorrect usage, and result in damage to the instrument. Improper usage could also cause severe personal injury, death, unpredictable results, instrument malfunction, and premature wear of components shortening the lifetime of the instrument. Such actions may void your warranty. Keep this user manual and any other safety and operating instructions provided with the instrument in a safe place that is accessible to all users for reference.

If you are concerned about the safe use of the instrument, please contact your authorized Miltenyi Biotec service provider or call Miltenyi Biotec Technical Support.

Content

EN	Important safety information	7
DE	Wichtige Sicherheitsinformationen	17
ES	Información importante de seguridad	27
FR	Informations de sécurité importantes	37
IT	Importanti istruzioni di sicurezza	47
1	Introduction to the MACSQuant® Instrument	57
2	Installation	59
2.1	Components included in the delivery	59
2.2	Complementary components	60
2.3	Unpacking the MACSQuant® Instrument	61
2.4	Installation of the MACSQuant® Instrument	63
2.4.1	Installing the fluid bottles	64
2.4.2	Installation of the Single Tube Rack	68
2.4.3	Installation of the MACS® MiniSampler	69
2.4.4	Installation of the MACSQuant® Column	72
2.5	Installation checklist	75
3	Switching on/off the MACSQuant® Instrument	77
3.1	Switching on the MACSQuant® Instrument for the first time	77
3.2	Switching off the MACSQuant® Instrument	78
3.2.1	Manual shutdown	78
3.2.2	Automatic shutdown options	78
4	The MACSQuantify™ Software	79

5	Calibration of the MACSQuant® Instrument hardware	81
5.1	Calibration of the sample uptake unit	82
5.1.1	Calibration of the robotic needle arm	82
5.1.2	Calibration of the Single Tube Rack	85
5.1.3	Calibration of the Chill 5 Rack	86
5.1.4	Calibration of the Chill 96 Rack	87
5.1.5	Calibration of the MACSPlex Filter Plate	88
5.1.6	Calibration of MACS® Reagent Rack 4	88
5.2	Calibration of the rack detection	90
5.2.1	Single Tube Rack and acquisition button detection	90
5.2.2	Rack detection	91
6	Hardware monitor	93
6.1	Fluidics	94
6.2	Sample uptake unit	95
6.3	Lasers and detectors	96
7	MACSQuant® Live Support	97
8	Maintenance	99
8.1	Decontaminating the MACSQuant® Instrument	99
8.2	Exchanging fluid bottles	100
8.3	Exchanging the MACSQuant® Column	100
8.4	Exchanging the Hydrophobic Air Filters	100
8.5	Cleaning the washing station	101
8.6	Cleaning the uptake needle	101
8.7	Long-term storage of the MACSQuant® Instrument	101

9	Troubleshooting	103
9.1	Column leakage	103
9.2	Syringe pump leakage	104
9.3	MACS® MiniSampler does not move properly	104
9.4	Air bubbles during measurement or failure of event acquisition	104
9.5	Excessive debris is present in acquisition	104
9.6	Touchscreen remains dark	105
10	Technical data and specifications	107
10.1	Technical data and specifications of the MACSQuant® Instrument	107
10.2	Technical data and specifications of the MACS® MiniSampler	112
11	Technical Support	113
12	Limited warranty	115
13	Glossary	117

Important safety information

WARNING

The MACSQuant® Instrument is designed for safe use when installed correctly, operated by trained personnel, operated in accordance with general safety practices and the instructions in this user manual.

Warnings and precautions

The guidelines in this section explain the potential risks associated with operating this instrument and provide important supplemental safety information to minimize the risks. Follow the instructions carefully to protect yourself, others, and the equipment from potential hazards and create a safe work environment. Use this instrument only as specified by the manufacturer to avoid damage of equipment and injury to personnel.

Always follow local working area safety instructions and laboratory policies as well as the standards for health, safety, and prevention of accidents. Contact the local authority governing electrical power supply, building constructions, maintenance, or safety for more information about the safe installation and operation of the instrument.

Hazard levels

Signal words are used to identify safety and property damage messages. The following signal words are used throughout this user manual.

WARNING

or **WARNING!** indicates a potentially hazardous situation, which, if not avoided, could result in death or serious injury.

CAUTION

or **CAUTION!** indicates a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Safety symbols

The following symbols are used to highlight conditions that could cause injury to personnel or damage to equipment.

Safety alert. Risk of danger. The documentation needs to be consulted in all cases where this safety symbol is used, in order to find out the nature of the potential hazard and any actions that have to be taken.

Risk of electric shock.

Strong magnetic field.

Persons carrying pacemakers or electronic medical implants must maintain distance.

Risk of crushing and shearing of bodily parts.

Hazardous optical radiation.

Hazardous laser radiation.

Biohazard. Risk of contamination if potentially dangerous biological material is used.

Protective conductor terminal. Symbol is attached on the inside of the instrument. Information for service personnel.

On (Power supply).

Off (Power supply).

Safety labels

Please notice the hazard points and safety symbols of the MACSQuant Instrument.


Figure 1: Hazard sources and the position of the safety symbols on the front of the MACSQuant Instrument.


Figure 2: Position of the safety symbols on the back of the MACSQuant Instrument.

The above mentioned safety labels and safety markings must be kept clean and legible. Periodically inspect the safety labels and safety markings and replace them if they are not legible or perceptible from a safe viewing distance. Contact Miltenyi Biotec for replacement labels.

General safety instructions

If the instrument is not working properly and/or the displayed instructions or messages advise you to contact Technical Support, safe operation of the instrument is no longer possible. Immediately switch off and unplug the instrument from the power outlet, and contact an authorized Miltenyi Biotec service provider or Miltenyi Biotec Technical Support.

Electricity and fire hazards

Potential risks arising from electrical devices include electric shock, short, and overheating. Electric shock may lead to severe personal injury or death. An electric short and overheating may lead to the spread of fire.

Electronic equipment might emit sparks, which could ignite combustible vapors or material resulting in explosion or fire. Do not use the instrument in areas classified as hazardous locations; for example, in oxygen-laden environments.

The housing of the instrument is designed to reduce the risks of electric shock, electric short, and spread of fire. Except for the front access covers, **do not remove or penetrate any cover**. Only authorized personnel must remove the other covers of the instrument. Never push a foreign object through an opening into the instrument. Do not use the instrument if

- it is opened or disassembled,
- it has been dropped or damaged,
- it has damaged or broken parts,
- it has a damaged power cable,
- an object has entered through the ventilation slots,
- a foreign object has been dropped into the instrument.

If flames or smoke appear, immediately switch off the power supply, unplug the instrument from the electrical outlet, and contact an authorized Miltenyi Biotec service provider or Miltenyi Biotec Technical Support.

The instrument is intended for indoor use only. **If fluids enter the instrument, it may lead to electric shorts, electric shock, or fire.** Take special care while handling fluids. Do not allow fluids to enter the interior of the instrument. Protect the instrument against accidental spillages and splashes. Clean up spillages immediately. Do not operate the instrument if liquids have spilled into the instrument. Do not use the instrument in a wet or damp location. Avoid areas with high humidity or condensation. Moving the instrument from a cold environment –

such as a cold room at 5 °C – to room temperature may cause condensation inside the instrument. In such cases, wait for the instrument to dehumidify before operating it. Unplug the instrument from the electrical outlet before cleaning. Do not use liquid or aerosol cleaning agents; always use a damp cloth.

Ensure adequate air circulation in the room where the instrument is operated. In the absence of adequate circulation, ambient air may not cool the instrument to acceptable operating temperatures. Allow sufficient air circulation around the instrument – at least 15 cm on all sides – during operation to ensure adequate cooling. Do not place the instrument near radiators, heat registers, stoves, or other equipment (including amplifiers) that generate heat. Protect the instrument from direct sunlight. Do not cover the slots and openings of the instrument as these are provided for ventilation, and protect the instrument from overheating. Do not place the instrument in a built-in rack or similar confined spaces unless the space has been specifically designed to provide proper ventilation. Follow the mounting instructions for the instrument.

Ensure that the main switch as well as the connector for the power cable are easily accessible and located as close to the operator of the instrument as possible. If it is necessary to disconnect the power supply, unplug the cable from the power outlet.

The instrument is equipped with a three-wire electrical grounding-type plug that has a third pin for grounding. This plug only fits into a grounded power outlet. This is a safety feature. Do not try to insert the plug into a non-grounded power outlet. If you cannot insert the plug into the outlet, contact your local electrician to replace the outlet.

The instrument should only be operated from a power source that meets the specifications mentioned on the product's electrical ratings label. If you have questions about the type of power source to use, contact your authorized Miltenyi Biotec service provider or local power company. Do not use extension cords or power strips. Do not overload an electrical outlet.

Only peripheral devices that comply to UL 60950 are allowed to be connected to the RS232 connector labeled COM 1. In addition, only original MACSQuant Instrument equipment should be attached to the connectors labeled External CAN. Only the MACSQuant Instrument fluid sensor cable should be attached to the Bottle Sensor connector. Only the 2D code reader that was delivered with the MACSQuant Instrument should be connected to the RS232/BCR connector. External laser devices connected to the connector labeled RS232/BCR have to comply with the standard IEC 60825-1. In order to ensure electromagnetic compliance connect peripheral devices only with double-shielded cables and only use the power cable and the keypad included in delivery. Only use connector cables less than 3 m in length. The keypad shall only be used with MACSQuant Instrument models.

Strong magnetic field

The MACSQuant Instrument is equipped with an extremely powerful magnet. **There is a risk of severe personal injury for persons carrying pacemakers, brain shunts, or electronic medical implants.** Keep any magnetic information carriers (such as credit cards, magnetic tapes, and floppy disks), any electronic equipment (such as hearing aids, pacemakers, brain shunts, measuring and control instruments, computers, and watches), and magnetizable tools and objects at a distance of at least 20 cm from the magnet cover. These items may be affected or damaged by the magnetic field.

Mechanical hazards

Moving and revolving parts are potential mechanical hazards. Do not open the front access covers while the instrument is in operation. Do not obstruct the movement of the robotic needle arm and accessory hardware during operation. There is a risk of puncturing bodily parts. Operate only with attached needle guard (see figure 2.10). **Keep fingers etc. away from all moving parts of the MACSQuant Instrument and accessory hardware, to avoid crushing or shearing injuries, or damage to the instrument.** Do not touch fluid pumps or adjust the tubing, while the instrument is in operation. Always switch off the instrument before adjusting any part of the fluidic system. Always stop or abort a procedure before handling accessory hardware, e.g., MACS® MiniSampler, or loading/removing tubes from the tube rack placed on the MiniSampler. **Do not circumvent any safety measures or devices.**

Optical radiation hazards

The MACSQuant Instrument is equipped with three Continuous-Wave Lasers (laser class 3B) for fluorochrome excitation. These lasers are secured by protective housing. Do not remove the protective housing. Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure. The MACSQuant Instrument is classified as a Class 1M laser product per standard IEC 60825-1:2007.

The instrument is equipped with four vertical cavity surface emitting lasers (VCSEL) for automated rack detection (Class 1M). The radiation is not visible. Do not view directly with optical instruments (e.g. lenses, magnifying glasses, and microscopes). Viewing the VCSEL port within 100 mm distance with optical instruments could be hazardous to the eye.

The instrument is also equipped with a 2D Code Reader which uses a visible semiconductor laser as a target pointer for adjusting the reading position, and with powerful light emitting diodes (LEDs) for illuminating the reading area and the fluid bottles as well as the waste bottle. Do not remove the bottle holders. According to

the international standard IEC 62471 this lamp system has an exposure hazard value (EHV) of 0.91 and is in excess of the Exempt Risk Group. The hazard distance (HD) for the Exempt Risk Group is 61 cm. The hazard distance for Risk Group 1 is 20 cm.

Do not look directly at laser or LED radiation or reflected laser or LED radiation from a mirrored surface. Otherwise, eye injury may result. Do not intentionally direct the laser beam at others. Be careful of the path of the laser beam or reflection from a mirrored surface. Take care during installation of the MACSQuant Instrument that the path of the laser beam is not at the same height as that of the human eye during operation. Do not allow water, oil, dust, or other foreign substances to stick to 2D Code Reader aperture window. This may cause read errors. Be sure to stop the laser emission before cleaning the scanner. Otherwise, exposure to the laser may cause eye injury. Use a soft, dry cloth to wipe any substances from the scanner. Do not use alcohol or other cleaning substance.

Do not disassemble, modify or remove the installed laser or LED radiation sources or their mounting brackets. The laser or LED radiation sources do not automatically stop emitting when disassembled. Radiation of disassembled units may lead to eye injuries.

Chemical and biological hazards

Depending on the biological material used, contamination or infection may lead to severe personal injury or death. All clinical samples must be considered potentially infectious. **If biohazardous material is or has been used, the operator must choose and wear personal safety equipment** as indicated in the warnings and precautions for the particular substance. The above safety precaution also accounts for any hazardous chemical, including toxic or corrosive chemicals, acidic, or radioactive substances, that may be present in the sample. Wear protective gloves and clothing, as well as safety glasses to prevent contact of the substance with skin and eyes. Also protect mouth and nose as aerosols might leak from the system. Defective or inadequate safety equipment is hazardous. The instrument must be operated in a safety hood if hazardous or unknown materials are being processed. If hazardous material has been used or spilled, care must be taken to thoroughly decontaminate the instrument. It is strictly prohibited to continue to handle contaminated accessories or parts of the instrument.

Tubes and any other consumables that were in contact with biohazardous material must be decontaminated before disposal. Liquid waste must be autoclaved or decontaminated using a disinfectant that is appropriate for the specific pathogen, e.g., 10% bleach, isopropyl alcohol, or 70% ethanol.

All liquid and solid waste must be considered biohazardous and, therefore, must be handled taking universal laboratory precautions. Waste disposal must be in accordance with any local regulations.

Always inspect the fluidics system (complete tubing set, reservoirs, bottles and their closures, valves, columns, diluters, peristaltic pumps, and uptake needle) before switching on the instrument. **If leakage has been detected, replace all damaged parts before switching on the instrument. If damaged parts cannot be replaced, unplug and do not use the instrument.** Failure of parts containing biohazardous material or liquids that have been in contact with such material could cause a hazard.

Servicing, transport, and disposal

Unless otherwise specifically noted in this user manual or other Miltenyi Biotec documentation, **do not service the instrument yourself.** Servicing and repair must be performed by qualified service personnel. Improper or incorrect unpredictable results, cause instrument malfunction or damage, as well as premature wear and reduced life of the instrument. It may also void your warranty.

When replacement or spare parts are required, ensure that the service provider uses only genuine Miltenyi Biotec parts or third-party parts specified and recommended by Miltenyi Biotec. Using unauthorized parts can cause malfunction of the instrument and impair results. Miltenyi Biotec does not honor any warranty or accept any responsibility for instrument failure or damages resulting from the use of inappropriate parts. After completing any service or repair, ensure that your authorized Miltenyi Biotec service provider performs all safety checks as required by the repair procedure to ensure that the instrument is operating correctly.

Only use options and upgrades recommended by Miltenyi Biotec.

Inquire with your local Miltenyi Biotec representative about Miltenyi Biotec's extensive instrument service and support arrangements, or refer to **www.miltenyibiotec.com/support**.

The instrument should be transported with care in packaging specified by Miltenyi Biotec. Internal damage can occur if the instrument is subjected to excessive vibration or dropped. If the instrument needs to be shipped back to the manufacturer for service, decontaminate the instrument to remove any hazardous material prior to shipment. If you have questions regarding proper decontamination or shipment, contact Miltenyi Biotec Technical Support for assistance.

Waste of Electrical and Electronic Equipment (WEEE) customer information

Please dispose of your end-of-life Miltenyi Biotec products in accordance with the applicable WEEE and hazardous waste disposal legislation, which may differ by country or region.

Electrical and electronic equipment may contain hazardous substances which may have a serious detrimental effect on the environment and/or human health. That is why all equipment must be specifically collected and treated by

designated waste facility centres and by qualified WEEE compliance schemes. By ensuring that you dispose of your unwanted electrical and electronic equipment according to the applicable WEEE and hazardous waste disposal legislation, you are helping to preserve our natural resources and protect human health.

Miltenyi Biotec is committed to protecting the environment. Miltenyi Biotec offers product end-of-life return programs in many countries, and partners with licensed WEEE compliance schemes throughout the world. Miltenyi Biotec lets you recycle your end-of-life Miltenyi Biotec equipment free of charge. The terms and availability of this offer vary by geography because of differences in regulatory requirements. Please note that, depending on the type and use of your equipment, additional requirements may apply.

For more information, or if you wish to dispose of your end-of-life Miltenyi Biotec equipment, please contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support.

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MACSQuant Instruments are compact flow cytometers that combine multisample and multiparameter sample analysis with unrivaled ease of use. With processing rates of up to 10,000 events per second, they facilitate also absolute quantitation of cell populations. MACS® Antibodies and Cell Separation Reagents are designed for easy and convenient handling with MACSQuant Instruments. Miltenyi Biotec is leading provider of comprehensive solutions for flow cytometry.


Figure 1.1: Front view of MACSQuant Instrument.

The key accessible components of the MACSQuant Instrument are:

- An integrated computer controlled by a touchscreen runs the MACSQuantify™ Software, which is used for sample acquisition as well as for data analyses.
- A fully automated computer-controlled robotic needle arm delivers samples to the sample injection port, which is located in the washing station. The uptake needle as well as the sample injection port are automatically washed during and after cell processing for your convenience.
- The MACS Single Tube Rack is used for the analysis of individual samples. It offers the start of your analysis by the push of a button.
- The optional MACS MiniSampler allows fully automated multiple sample processing. It holds different tube racks that are automatically recognized by a 2D code reader, as well as the MACS Reagent Rack 4.
- A color-coded and easy to access fluidics system supplies the MACSQuant Instrument with all buffers and solutions required during and after operation. The sheath particle filter prevents any debris from entering the fluidics system.
- Bottles for Running Buffer, Washing Solution, and Storage Solution are placed in bottle holders. Liquid waste is collected in a separate waste bottle. All fluid levels are constantly monitored and visualized by the illumination of the fluid bottles.
- The MACS Enrichment Unit with the optional MACSQuant Column allows easy analysis of particularly challenging cell samples, such as rare cell populations in high volumes of buffer, or precious low volume samples.

2

Installation

Read the instructions in the chapter **Important safety information** and ensure that your site is properly prepared before continuing with installation and assembly. Incorrect installation can lead to the spread of fire, explosion, the risk of electric shock, or biohazards.

The following section describes how the MACSQuant Instrument is unpacked and installed for first use. Please note that a qualified Miltenyi Biotec trained representative should perform the unpacking and installation of the MACSQuant Instrument.

2.1 Components included in the delivery

Component	Position in transportation box
MACSQuant Instrument	Main storage compartment
Fluid sensor cable module	Top storage compartment
Power cable	Top storage compartment
Fluidics tubing (6x)	Top storage compartment
Bottle closures with sensors (4x)	Top storage compartment
Single Tube Rack	Top storage compartment
User manual	Top storage compartment
Empty bottles with closures	In bottle holders
Keyboard	Front compartment
Webcam	Top storage compartment
Hydrophobic Air Filters	Top storage compartment
MACSQuant Calibration Beads	Top storage compartment
MACSQuant Starting Buffer Kit	Extra packaging

2.2 Complementary components

Complementary components can be ordered separately.

Component	Description	Order no.
MACS MiniSampler	Together with the MACS Reagent Rack 4 for automated multisample processing.	130-092-535
Chill 5, 15, 50 Rack Set	Set of cooling racks	130-097-038
Chill 5 Rack	Cooling rack for 24×5 mL tubes	130-092-951
Chill 15 Rack	Cooling rack for 15×15 mL tubes	130-092-952
Chill 50 Rack	Cooling rack for 6×50 mL tubes	130-092-953
Chill 96 Rack	Cooling rack for 96-well plates	130-094-459
MACSQuant Columns	For use with the MACS Enrichment Unit	130-094-458
Hydrophobic Air Filter	To vent fluid bottles	130-090-385
MACS Laminar Hood Plate	Flat working platform for the MACSQuant Instrument	130-093-246
MACSQuant Remote Port	For connection to liquid handling systems	130-097-215
Buffer Supply System	For use of larger fluid containers	130-098-628
MACSQuant Calibration Beads	For instrument calibration	130-093-607
MACSQuant Running Buffer	6×1.5 L	130-092-747
MACSQuant Washing Solution	6×1.5 L	130-092-749
MACSQuant Storage Solution	6×1.5 L	130-092-748
MACSQuant Starting Buffer Kit	4×1.5 L Running Buffer 1×1.5 L Washing Solution 1×1.5 L Storage Solution	130-094-190
MACSQuant Washing & Storage Solution Kit	3×1.5 L Washing Solution 3×1.5 L Storage Solution	130-092-801
MACS Bleach Solution	6×1.0 L	130-093-663
MACSQuant IQ/OQ	Installation and operational qualification service for MACSQuant Instruments	160-001-711
Full Service	Service contract for the MACSQuant Analyzer or MACSQuant Analyzer 10	160-001-287
Full Service VYB	Service contract for the MACSQuant VYB	160-001-637

2.3 Unpacking the MACSQuant® Instrument

WARNING! If there is any damage, do not use the instrument but contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support. The operation of a damaged instrument may lead to the exposure of mechanical hazards, optical radiation hazards, electric shock, or the spread of fire. Read the chapter **Important safety information** before continuing with installation and assembly. Read through the following instructions carefully before commencing the installation procedure. Before opening the transportation box, check for any visible external damage to the box. Check also to see if the shock and position indicators (if present) suggest incorrect transportation of the instrument.

The MACSQuant® Instrument is designed to fit perfectly onto benchtops, into standard sized laminar flows, into safety cabinets, or in connection to fully automated liquid handling systems. It should be installed on a stable, flat and vibration-free surface. The operating environment should be dust-free, sufficiently ventilated, and free from sources of electromagnetic radiation. In order to ensure a flat surface in laminar flow hoods, the MACSQuant Instrument can be placed on a MACS Laminar Hood Plate.

- 1** Open the transportation box. The top storage compartment holds the MACSQuant Instrument user manual, cardboard boxes and bags, which contain components included in the delivery as well as any complementary components.
- 2** Take out all boxes and bags from the top storage compartment and remove the foam.
- 3** Take out the keyboard from the front protection foam. Remove the front protection foam and the two side protection foams.

- 4 Lift the MACSQuant Instrument by holding it at the base of the orange bottle holders at both sides of the device. **CAUTION!** Two persons together should lift the MACSQuant Instrument. Take care to avoid personal strain or injury. The instrument is heavier at the front and should be stabilized while lifting it. Miltenyi Biotec accepts no liability for potential injuries sustained during lifting and/or movement of the device.

Figure 2.1: Two persons are required to lift the MACSQuant Instrument.

- 5 Transfer the instrument onto a stable working platform and remove the plastic bag surrounding the device. **WARNING!** Take into consideration that the instrument requires adequate air circulation for heat exchange and cooling. Read the chapter **Important safety information** to avoid the risk of overheating. Do not place the instrument on an unstable table, cart, stand, tripod, or bracket. As a consequence, the instrument might fall down, which may cause serious bodily harm and/or serious damage to the instrument. Use only a table, cart, stand, tripod, or bracket that can easily support a weight of >50 kg.
- 6 Carefully remove the uptake needle from the foam packaging.

Figure 2.2: The uptake needle and the stylus are stored in the foam packaging.

- 7 Place the uptake needle into its guiding at the robotic needle arm. Ensure that the tubing connected to the uptake needle can move freely when the robotic needle arm extends, or when the uptake needle moves into the sample uptake position. Remove the protective cover of the needle.

Figure 2.3: The uptake needle is placed into its guiding at the robotic needle arm.

- 8 Remove stylus from packaging and place in holder attached to touchscreen.
- 9 Adjust the angle of the touchscreen in an upright position in order to remove the foam packaging.
- 10 Proceed to chapter 2.4 **Installation of the MACSQuant Instrument.**

2.4 Installation of the MACSQuant® Instrument

WARNING! Before operating the MACSQuant® Instrument for the first time, carefully read the user manual and contact your local Miltenyi Biotec representative for assistance. When delivered, the MACSQuant Instrument fluidics system is delivered dry without Storage Solution. When working with biohazardous samples, it is recommended to fill the waste bottle with 100 mL of disinfectant, e.g. MACS Bleach Solution, before use. Empty bottles from the Running Buffer, Washing Solution, or Storage Solution may be used as waste bottles only with the appropriate safety label on it. For proper disposal, please follow local regulations and carefully read the chapter **Important safety information.**

2.4.1 Installing the fluid bottles

Operating the MACSQuant Instrument requires three fluid bottles for Running Buffer, Washing Solution, and Storage Solution, as well as a waste bottle. Always operate the instrument with ready-to-use MACS Buffers and solutions only. The bottles, closures, tubing, and fluid sensor cables are color coded for easier handling.

Bottle	Symbol
Running Buffer (blue)	
Washing Solution (green)	
Storage Solution (black)	
Waste (red)	

Table 2.1: Symbols and color-coding of fluid bottles.

- 1 Install one fluid bottle at a time. Place a new bottle into the orange bottle holder. Please note the corresponding color-coding.
- 2 Unscrew the lid of the fluid bottle and replace it with the appropriate bottle closure delivered with the MACSQuant Instrument.
WARNING! Do not unscrew the bottle lid until the bottle is placed in the bottle holder to avoid spilling on the instrument and risk of electric shock.

- 3 Remove the caps from the bottle distribution block and connect the tubing to the appropriate color-coded fluid port of the bottle distribution block on the back of the MACSQuant Instrument.

Figure 2.4: Connect the color-coded tubing to the according position of the bottle distribution block.

- 4 Note the color coding and connect the other end of the tubing to the respective bottle closure.

Figure 2.5: Connect the according color-coded tubing and sensor cables to the appropriate connectors on the waste bottle closure and the Storage Solution bottle closure.

- 5 Note the color coding and connect each fluid sensor cable to the respective bottle closure.

Figure 2.6: Connect the according color-coded tubing and sensor cables to the appropriate connectors on the Running Buffer bottle closure and the Washing Solution bottle closure.

- 6 Connect the fluid sensor cable to the port labeled **Bottle Sensor** at the back of the MACSQuant Instrument and fasten securely.

Figure 2.7: Connect the sensor cables to the sensor cable port on the back of the MACSQuant Instrument.

- 7 Attach the cable coming from the 2D code reader to the port labeled **RS232/BCR** at the back of the MACSQuant Instrument and fasten securely.

Figure 2.8: Connect the cable of the 2D code reader to the RS232/BCR port of the MACSQuant Instrument.

- 8 Connect Hydrophobic Air Filters to the appropriate connectors on the bottle closures to keep buffers sterile. Avoid any contact of Hydrophobic Air Filters with fluids as this may cause clogging of the filters.
- 9 Connect a Hydrophobic Air Filter to the appropriate connector of the fluidics system on the front of the MACSQuant Instrument.
- 10 Connect a Hydrophobic Air Filter to the appropriate connector behind the sheath particle filter at the back of the MACSQuant Instrument. Please note that the white lid on the top of the sheath particle filter may need to be unscrewed in order to purge the system.

Figure 2.9: Purge the system by opening the white lid on top of the sheath particle filter.

- 11 Check if the needle guard is attached to the robotic needle arm and attach it if needed.

Figure 2.10: The needle guard is attached to the robotic needle arm in order to protect against the risk of puncturing bodily parts.

2.4.2 Installation of the Single Tube Rack

The MACSQuant Instrument is delivered with the Single Tube Rack by default. The Single Tube Rack has been designed for compatibility with all standard flow cytometry tubes (1.5 mL, 2 mL, and 5 mL) and offers the start of your analysis by the push of a button. Gently insert the Single Tube Rack into the corresponding slots located at the front of the MACSQuant Instrument. The rack should click into place. The included CAN plug must be attached to the back of the MACSQuant Instrument in the port labeled **External CAN** (see figure 2.15).

In order to remove the Single Tube Rack, gently pull it by holding it between two fingers. Please note that the orange button will start a new measurement of the MACSQuant Instrument and does not unlock the Single Tube Rack.

Figure 2.11: Install or remove the Single Tube Rack by simply plugging or unplugging.

2.4.3 Installation of the MACS® MiniSampler

The MACSQuant Instrument is optionally delivered with the MACS® MiniSampler, the MACS Reagent Rack 4, and different cooling racks, which can be pre-cooled for 3 to 4 hours at 2–8 °C. Please note that the cooling racks shall not be pre-cooled in a freezer. Once installed, the MACS MiniSampler is automatically recognized by the MACSQuant Instrument. Each tube rack has a barcode on the rear side that is detected upon placement of the rack on the MACS MiniSampler. Before the MACS MiniSampler can be installed, please remove the Single Tube Rack.

- 1 Remove the transparent protection foil from the lens of the rack detector. Note the positions of the MACS MiniSampler slot located at the front of the MACSQuant Instrument and the MACS MiniSampler guiding.

Figure 2.12: MACS MiniSampler slot at the front of the MACSQuant Instrument and the corresponding guiding on the MACS MiniSampler.

- 2 Tilt the MACS MiniSampler and slide the guiding into the receiving slot until resistance is met; lower the MACS MiniSampler to a horizontal position in order to lock the MACS MiniSampler.


Figure 2.13: Installation of the MACS MiniSampler.

- 3 Open the front door of the MACSQuant Instrument and fix the MACS MiniSampler with the according screw.

Figure 2.14: Fixing the MACS MiniSampler with the according screw.

- 4 Close the front door of the MACSQuant Instrument and ensure that the MACS MiniSampler is completely inserted and secured.
- 5 Place the MACS MiniSampler cable underneath the MACSQuant Instrument and connect it to the port labeled **External CAN** at the back of the instrument.

Figure 2.15: The MACS MiniSampler cable is attached to External CAN port.

- 6 Install the MACS Reagent Rack 4 onto the MACS MiniSampler into the left recess by snapping the engagement hook into the undercut.


Figure 2.16: The MACS Reagent Rack 4 is installed on the MACS MiniSampler.

- 7 Place a MACS Chill Rack onto the MACS MiniSampler into the right recess ensuring that the rack barcode is facing the MACSQuant Instrument. If automatic recognition of the tube rack fails, the instrument will display a screen for manual selection of the tube rack. Before confirming the choice, ensure that the rack is placed correctly into the recess.

Figure 2.17: The Chill 5 Rack is installed on the MACS MiniSampler.

- 8 Note the position of the lid guiding at both sides of the MACS MiniSampler and attach the protective cover. Keep the protective cover of the MACS MiniSampler attached to the instrument and closed when operating.

Figure 2.18: The protective cover is installed on the MACS MiniSampler.

2.4.4 Installation of the MACSQuant® Column

WARNING! Always wear protective gloves and eyewear in order to protect against potential biohazard exposure.

The MACSQuant® Instrument is delivered with a column substitute, which can be exchanged for a MACSQuant Column. MACSQuant Columns are specifically designed for the use with the MACS Enrichment Unit, which is integrated into the MACSQuant Instrument. Magnetically and fluorochrome labeled cells are pre-enriched in the MACSQuant Column prior to cell analysis allowing flow cytometric analysis of rare cell populations in a fraction of the time. Up to 5×10^6 magnetically labeled cells can be enriched and analyzed in a single step. The entire process of sample labeling, pre-enrichment and flow cytometric analysis is fully automated by the MACSQuantify Software. The MACSQuant Column can be used for up to three months before it needs to be replaced. Please keep the column substitute; it can be reinstalled if the MACSQuant Column is no longer needed.

- 1 Open the front access cover of the MACSQuant Instrument in order to reach the MACS Enrichment Unit and the column substitute/MACSQuant Column.

- Using both hands, hold the top and bottom of the column substitute and pull gently but firmly to remove it from its slot in the MACS Enrichment Unit.

Figure 2.19: Removing the column substitute from the MACS Enrichment Unit.

- Hold the column in one hand and gently unscrew the upper column connector counter-clockwise. Tilt the column downwards to empty any fluid into a paper towel placed underneath.
- Unscrew the bottom column connector to remove the column.
- Screw the bottom column connector into the bottom end of the new MACSQuant Column. Please note the right orientation of the column. The top end of the MACSQuant Column has a clearly visible inlet, while the metal balls reach further to the bottom end of the column. The column should be inserted only with the visible inlet at the top position.

Figure 2.20: Installation of the bottom connector to the MACSQuant Column.

- 6 Point the column towards the top of the device and screw in the top column connector.

Figure 2.21: Installation of the top connector to the MACSQuant Column.

- 7 Align the column so that the top column connector sits on the guide of the magnet cover. Push the column into the slot. A clear click can be perceived. Verify that the column is placed in the center of the magnet cover.
- 8 Close the front access cover of the MACSQuant Instrument.
- 9 Prime the new MACSQuant Column by selecting the option **Prime column** from the **Tools** tab in the **Sidebar** of the MACSQuantify Software. Please see also the MACSQuant Instrument software guide for more information.

To replace a MACSQuant Column follow the same procedure.

2.5 Installation checklist

Please ensure the correct installation of the MACSQuant Instrument by the following checklist:

- 1** Note the correct positioning of each fluid bottle and the waste bottle, recognizable by the color code and the symbols. This is crucial for successful analyses using the MACSQuant Instrument.
- 2** Ensure that the correct tubing and fluid sensor cable is attached to the corresponding bottle, and that all tubing connections as well as the bottle closures are fastened. If necessary, tighten loose connections using a wrench.
- 3** Ensure that the fluid bottles are filled and installed correctly. For safe operation of the MACSQuant Instrument, all fluid bottles must contain at least 150 mL of the respective solution. In order to prevent contamination of the fluidics system the use of non-sterile buffers and solutions is not recommended. Only use buffers and solutions supplied by Miltenyi Biotec for operation of the MACSQuant Instrument. A reproducible and optimal performance of the MACSQuant Instrument cannot be guaranteed when the instrument is operated using self-made buffers and/or solutions procured from another manufacturer.
- 4** Make sure that the waste bottle is empty and installed correctly.
WARNING! When working with biohazardous samples, it is recommended to fill the waste bottle with 100 mL of disinfectant (e.g. MACS Bleach Solution) before use. Empty bottles from the Running Buffer, Washing Solution, or Storage Solution may be used as waste bottles only with the appropriate safety label on it. For proper disposal, please follow local regulations and carefully read the chapter **Important safety information**.
- 5** Open the front door and check that the column substitute or the MACSQuant Column is installed correctly. Ensure that the tubes are securely fastened to the column and that no part of the visible tubing is pinched or obstructed.
- 6** Ensure that the needle guard is attached to the MACSQuant Instrument.
WARNING! There is a hazard of crushing, shearing, or puncturing bodily parts. Operate only with attached needle guard.
- 7** When all points of this installation checklist have been fulfilled, close all covers. Proceed to chapter 3 Switching on/off the MACSQuant Instrument.

3

Switching on/off the MACSQuant® Instrument

3.1 Switching on the MACSQuant® Instrument for the first time

- 1 Connect the power cord to the power socket on the rear panel of the MACSQuant® Instrument. The main power switch is located on the right side of the instrument in front of the bottle holders (I indicates on, O indicates off). Ensure that the main switch is in position O before connecting the power cord.
- 2 Plug in the power cord to a grounded power outlet. **WARNING!** The device is equipped with a three-wire electrical grounding-type plug. This is a safety feature to reduce the risk of electric shock. Do not insert the plug into a non-grounded power outlet. If you cannot insert the plug into the outlet or if you are not sure whether the outlet is grounded, contact your local electrician. The device should only be operated from a power source indicated on the product's electrical ratings label. If you have questions about the type of power source to use, contact your authorized Miltenyi Biotec service provider or local power company. Ensure that the main switch as well as the connector for the power cable are easily accessible and located as close to the operator of the instrument as possible. In case of danger disconnect the device from the main supply.
- 3 Switch on the MACSQuant Instrument to turn it into stand-by mode. The MACSQuantify Software will automatically power up if you touch the touchscreen.

3.2 Switching off the MACSQuant® Instrument

Please note that before switching off the MACSQuant® Instrument shall be shut down. Refer to the MACSQuant Instrument software guide for an introduction to the MACSQuant Instrument user interface.

3.2.1 Manual shutdown

If no further measurements are required, shut down the instrument as follows:

- 1 Click the **Main instrument control** button from the **Top menu** of the MACSQuantify Software and select instrument off. This automatically starts a washing protocol that lasts for approximately seven minutes. This protocol includes an incubation of the Washing Solution in the fluidics followed by the flushing of the Washing Solution and replacement with the Storage Solution. With the **Main instrument control** button the MACSQuant Instrument can also be switched to data analysis mode, which initiates the washing protocol, too, but enables data analysis using the software interface.
- 2 After the end of the washing protocol the MACSQuant Instrument is in stand-by mode. It can be securely switched off now or, if desired, rebooted by touching the touchscreen.

3.2.2 Automatic shutdown options

An automatic shutdown procedure is set as default, which prevents the non-essential illumination of lasers that can shorten diode lifespan. Modification of the default settings of the automated shutdown options is possible under **Edit > Options (default) > Software > Timers** from the **Menu bar** of the MACSQuantify Software for administrators only. The **Standby timer** is the amount of idle time (no fluidic operation) before automated shutdown begins. The **Shutdown timer** is the incubation time during the shutdown procedure, which is recommended to be at least 5 minutes. The **Needle priming time** ensures that air is not introduced during measurement due to evaporation while the uptake needle is idle. The recommended priming time is 20 seconds. For details on adjusting the automated shutdown options, please refer to the MACSQuant Instrument software guide and/or the online help of the MACSQuant Instrument.

4

The MACSQuantify™ Software

The MACSQuantify™ Software is installed on every MACSQuant Instrument and can be used both for sample acquisition as well as for data analyses. The integrated computer offers various ports to connect external devices (e.g., external keyboard and mouse, USB flash device, webcam, external display, ethernet network cable). The MACSQuantify Software can also be installed onto personal computers (PCs) to allow data analysis independent from the MACSQuant Instrument or from other flow cytometers. The MACSQuantify Software runs on Windows XP® (SP3) or higher operating system with a minimum of 1 GB random access memory (RAM) only. The MACSQuantify Software when installed onto a PC will operate for a trial period of 30 days without registration. Please contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support to register. For updates of the MACSQuantify Software on the MACSQuant Instrument, if available, follow the instructions for software updates as directed by the Miltenyi Biotec representative or Miltenyi Biotec Technical Support.

For detailed instructions for use of the MACSQuantify Software, please refer to the MACSQuant Instrument software guide and/or the online help of the MACSQuant Instrument.

5

Calibration of the MACSQuant® Instrument hardware

Read the instructions in the chapter **Important safety information** before operation of the instrument.

When using the MACSQuant® Instrument for the first time, it is necessary to calibrate the hardware. For optimal operation, it is necessary to calibrate also the photomultiplier tubes (PMTs) of the instrument. It is highly recommended that hardware calibration is only performed by an administrator who has been trained by Miltenyi Biotec. In order to calibrate the MACSQuant Instrument the user must be familiar with the MACSQuantify Software. Refer to the MACSQuant Instrument software guide and/or the online help of the MACSQuant Instrument for an introduction to the MACSQuant Instrument user interface and a detailed description of the instrument settings.

5.1 Calibration of the sample uptake unit

Click on the **Tools** tab in the **Sidebar** of the MACSQuantify Software and check the **Calibrate uptake unit** box. Red closed circles indicate that the uptake components are not calibrated.


Figure 5.1: Selection of Tools tab for uptake unit calibration.

5.1.1 Calibration of the robotic needle arm

The robotic needle arm moves between the samples and the washing station with its sample injection port and washing port.

- 1 Ensure that the robotic needle arm can freely move and that no object is obstructing it.

- 2 Click **Start calibration** under the heading **Needle arm calibration relative to washing station**. The robotic needle arm will move toward the washing station before being inserted into the sample injection port. When the uptake needle pops out of the robotic needle arm holder, a system failure is caused and a prompt screen to reinitialize appears. Reinsert the uptake needle, adjust the robotic needle arm appropriately, and confirm the prompt screen by pressing **OK**.


Figure 5.2: Calibration of the uptake unit is underway. The yellow closed circle indicates that the robotic needle arm may be freely moved for manual calibration.

- 3 Manually adjust the robotic needle arm in order to position the uptake needle over the center of the sample injection port.

- 4 Gently lower the robotic needle arm into the sample injection port until the uptake needle makes first contact with the bottom of the orifice. The uptake needle should always remain perpendicular to the horizontal plane, i.e., should not be arched or inserted diagonally.


Figure 5.3: Inserting the robotic needle arm into the sample injection port of the washing station.

- 5 When the uptake needle stands on the bottom of the sample injection port, lower the uptake needle holder by another extra 5 mm.

Figure 5.4: Lower the robotic needle arm holder by another 5 mm.

- Click **Save calibration**. The closed green circle indicates that the calibration is completed. Click **Test** to confirm that the correct coordinates are saved.


Figure 5.5: Coordinates of the robotic needle arm are saved.

5.1.2 Calibration of the Single Tube Rack

- Attach the Single Tube Rack with desired sample tube to the MACSQuant Instrument as described in chapter 2.4.2
Installation of the Single Tube Rack.
- Click the **Experiment** tab in the **Sidebar** of the MACSQuantify Software and select the **Single Tube Rack** format from the **Rack** pull-down menu. If the adjacent checkbox is activated the rack will be detected automatically. If rack detection has not previously been calibrated, deactivate this checkbox before proceeding.


Figure 5.6: Selecting Single Tube Rack for rack calibration.

- 3 Click **Start calibration**. The robotic needle arm will automatically move forward and insert the uptake needle into the single tube. The uptake needle should be positioned at the center of the tube (i.e. equidistant from the tube edges) almost touching the tube bottom (i.e. only a fraction of a millimeter from the bottom of the tube). To check the uptake needle position, gently wiggle the tube to ensure that there is a small amount of movement.
- 4 If this is not the case, carefully adjust the robotic needle arm accordingly.
- 5 Click **Save calibration**. The closed green circle indicates that the calibration is completed. Click **Test** to confirm that the correct coordinates are saved.

5.1.3 Calibration of the Chill 5 Rack

For the correct usage of the MACS MiniSampler, please calibrate the Chill 5 Rack. If the Chill 5 Rack is calibrated correctly, the Chill 15 Rack and the Chill 50 Rack are also automatically calibrated.

- 1 Attach the MACS MiniSampler to the MACSQuant Instrument as described in chapter 2.4.3 **Installation of the MACS MiniSampler**.
- 2 Click the **Experiment tab** in the **Sidebar** of the MACSQuantify Software and select the **Chill 5 tube rack** format from the **Rack** pull-down menu. If the Chill 5 Rack is not selected, the MACSQuant Instrument will display a prompt screen **Specified rack must be specified in the experiment settings**.


Figure 5.7: Selecting Chill 5 Rack using the Experiment tab.

- 3 Place a Chill 5 Rack onto the MACS MiniSampler and load an empty 5 mL tube in the **position D6**.
- 4 Click **Start calibration**. The robotic needle arm will automatically move forward and insert the uptake needle into the tube in the **position D6** of the Chill 5 Rack. The uptake needle should be positioned at the center of the tube (i.e. equidistant from the tube edges) almost touching the tube bottom (i.e. only a fraction of a millimeter from the bottom of the tube). To check the uptake needle position, gently wiggle the tube to ensure that there is a small amount of movement.

Figure 5.8: Calibration of the Chill 5 Rack by defining the position D6.

- 5 If this is not the case, carefully adjust the robotic needle arm accordingly.
- 6 Click **Save calibration**. The closed green circle indicates that the calibration is completed. Click **Test** to confirm that the correct coordinates are saved. The MACSQuant Instrument will automatically test selected positions.

5.1.4 Calibration of the Chill 96 Rack

- 1 Click the **Experiment tab** in the **Sidebar** of the MACSQuantify Software and select the **96-well** format from the **Rack** pull-down menu. If the 96-well is not selected, the MACSQuant Instrument will display a prompt screen **Specified rack must be specified in the experiment settings**.


Figure 5.9: Selecting a 96-well rack for rack calibration.

- 2 Place a Chill 96 Rack onto the MACS MiniSampler and load an empty 96-well plate.

- 3 Click **Start calibration**. The robotic needle arm will automatically move forward and insert the uptake needle into the rack **position H12** of the Chill 96 Rack. The uptake needle should be positioned at the center of the well (i.e. equidistant from the well edges) almost touching the well bottom (i.e. only a fraction of a millimeter from the bottom of the well). To check the uptake needle position, gently wiggle the plate to ensure that there is a small amount of movement.
- 4 If this is not the case, carefully adjust the robotic needle arm accordingly.
- 5 Click **Save calibration**. The closed green circle indicates that the calibration is completed. Click **Test** to confirm that the correct coordinates are saved. The MACSQuant Instrument will automatically test selected positions.

5.1.5 Calibration of the MACSPlex Filter Plate

- 1 Click the **Experiment tab** in the **Sidebar** of the MACSQuantify Software and select the **MACSPlex Filter Plate** format from the **Rack** pull-down menu. If the MACSPlex Filter Plate is not selected, the MACSQuant Instrument will display a prompt screen **Specified rack must be specified in the experiment settings**.
- 2 Place a Chill 96 Rack onto the MACS MiniSampler and load an empty MACSPlex Filter Plate.
- 3 Click **Start calibration**. The robotic needle arm will automatically move forward and insert the uptake needle into the rack **position H12** of the MACSPlex Filter Plate. The uptake needle should be positioned at the center of the well (i.e. equidistant from the well edges) almost touching the well bottom (i.e. only a fraction of a millimeter from the bottom of the well). To check the uptake needle position, gently wiggle the plate to ensure that there is a small amount of movement.
- 4 If this is not the case, carefully adjust the robotic needle arm accordingly.
- 5 Click **Save calibration**. The closed green circle indicates that the calibration is completed. Click **Test** to confirm that the correct coordinates are saved. The MACSQuant Instrument will automatically test selected positions.

5.1.6 Calibration of MACS® Reagent Rack 4

- 1 Place the MACS® Reagent Rack 4 onto the MACS MiniSampler as described in chapter 2.4.3 **Installation of the MACS MiniSampler**.
- 2 Load an empty reagent vial in **position 1** of the MACS Reagent Rack 4 and close the protective cover.
- 3 Click **Start calibration**. A dialog box will prompt you to ensure that the MACS Reagent Rack 4 is indeed correctly placed and that the Single Tube Rack is removed.


Figure 5.10: Dialog box prompting you to ensure that the MACS Reagent Rack 4 is correctly installed.

- 4 Confirm the dialog box by clicking **OK**. The robotic needle arm will automatically move forward and insert the uptake needle into the reagent vial **position 1**. If this is not the case, carefully adjust the robotic needle arm accordingly.

Figure 5.11: Calibration of the Reagent Rack 4.

- 5 Click **Save calibration**. The closed green circle indicates that the calibration is completed. Click **Test** to confirm that the correct coordinates are saved. The MACSQuant Instrument will automatically test selected vial positions.

5.2 Calibration of the rack detection

5.2.1 Single Tube Rack and acquisition button detection

- 1 Click on the **Tools** tab in the **Sidebar** of the MACSQuantify Software and check the **Calibrate rack detection** box. Red closed circles indicate that the components are not calibrated.
- 2 Click **Start calibration** under the heading **Single Tube Rack and acquisition button detection**.
- 3 When prompted, remove the Single Tube Rack by gently pulling the component away from its docking-adaptor as described in chapter 2.4.2 **Installation of the Single Tube Rack**. Please note that the orange button is used only for acquisition; it does not release the Single Tube Rack from the MACSQuant Instrument.
- 4 When prompted, reattach the Single Tube Rack to the MACSQuant Instrument.
- 5 When prompted, push the orange sample acquisition button of the Single Tube Rack until the calibration is completed. The closed green circle indicates that the calibration is completed.


Figure 5.12: Successful completion of single rack detection.

5.2.2 Rack detection

This step is required for the automatic detection of Chill 5, 15, and 50 Racks by the MACSQuant Instrument.

- 1 Attach the MACS MiniSampler to the MACSQuant Instrument as described in chapter 2.4.3 **Installation of the MACS MiniSampler**.
- 2 Click the **Experiment** tab in the **Sidebar** of the MACSQuantify Software and select the **Chill 5 tube rack** format from the **Rack** pull-down menu. If the adjacent checkbox is activated the rack will be detected automatically. Deactivate this checkbox before proceeding.


Figure 5.13: Selecting Chill 5 Rack using the Experiment tab.

- 3 Place a Chill 5 Rack onto the MACS MiniSampler and click **Start calibration** under the heading **Rack detection**. The MACSQuant Instrument will move the Chill 5 Rack in order to check the 2D code reader for operability. The closed green circle indicates that the calibration is completed. After all calibrations are completed, the instrument can be switched into acquisition mode.


Figure 5.14: Rack detection is successfully completed.

6

Hardware monitor

Read the chapter **Important safety information** before operating this instrument. Never run the instrument if any errors are indicated within the Hardware monitor.

For detailed instructions for use of the MACSQuantify Software, please refer to the MACSQuant Instrument software guide and/or the online help of the MACSQuant Instrument. The hardware monitor can be accessed under the **View** menu. With the hardware monitor the real time status of the MACSQuant Instrument can be assessed, for example, to provide additional information in the case of error messages appearing on the screen. Here you find also the serial number of the MACSQuant Instrument and informations on the photomultiplier tubes (PMT).

6.1 Fluidics

The fluidics tab in the hardware monitor displays the pumps and valves in the live status, as well as the status of the fluid bottles.

Component	Further information
Fluid bottle monitor	Displays buffer/solution levels in real time
Pump monitor	Displays the status of the waste (W), air (A), and fill (F) pumps
Separation unit monitor	Displays the status of the MACSQuant Column and MACS Enrichment Unit
Valves 1–3 monitor	Displays the position of the valves for the MACSQuant Column: C = closed, O = open, green = in use
Valves 4–6 monitor	Displays the position of the general fluidics system valves
Syringe drive monitor	Displays the position of the dilutor
Sensor monitor	Displays the general system pressure and fluid reservoir levels.

Table 6.1: The different panels of the fluidics tab

The waste (**W**), air (**A**), and fill pump symbols (**F**) are illuminated in green when active. Active separator valves (**Valves 1-3**) are indicated in green, defective rotary valves in red. Valve status is indicated by o for open and c for closed. The fluid bottles in the hardware monitor as well as on the instrument will be illuminated in red when empty; the waste bottle will be illuminated in red when full.


Figure 6.1: Real-time Hardware monitor of the fluidic components.

6.2 Sample uptake unit

The sample uptake tab indicates the status of the robotic needle arm as well as whether the MACS MiniSampler is connected or not. The relative position of the robotic needle arm is listed in the lower left box as well as whether the Single Tube Rack is connected or not.


Figure 6.2: Real-time hardware monitor of the sample uptake arm.

6.3 Lasers and detectors

The lasers and detectors tab displays the status of the optical bench. It is possible to monitor the status of each laser fluorescence channel. The temperature, fan speed, PMT voltage, and annotated path of each laser is shown. A status overview of the optical bench is schematically represented. The bench temperature is kept between 33 °C and 37 °C. Therefore the fan speed is regulated automatically depending on the ambient room temperature and internal temperature of the MACSQuant® Instrument. Please note that the temperature for the lasers should be between 10 °C and 45 °C. In case of errors, please contact Miltenyi Biotec Technical Support.


Figure 6.3: Real-time Hardware monitor of the optical bench.

7

MACSQuant® Live Support

MACSQuant® Live Support is a real-time diagnostic service provided by Miltenyi Biotec Technical Support. Highly trained MACSQuant Specialists can be reached in real-time to assist with any queries you may have. A webcam is provided for communication with the MACSQuant Live Support, which is attached via USB connection to the MACSQuant Instrument. The webcam is installed to a free USB port at the back of the instrument and detected automatically by the software. During normal operation the webcam, however, should not be installed on the instrument. Please note that the MACSQuant Instrument must have network access to the internet in order to use the MACSQuant Live Support.

- 1 Select **MACSQuant Live Support** from the **Tools** tab in the **Sidebar** of the MACSQuantify Software to access remote assistance.
- 2 Complete all fields and detail any queries you may have using the **Message/Question** box.
- 3 Click **Submit** to commence MACSQuant Live Support.

HTML Help

Welcome to MACSQuant Live Support

To talk with a representative, please answer the following:

First name: (Required)

Last name: (Required)

Country / City: (Required)

Email address: (Required)

Phone number: (Required)

Message / Question:

Figure 7.1: Dialog box of the MACSQuant Live Support.

8

Maintenance

Read the chapter **Important safety information**. Do not service the equipment yourself. Service shall be carried out only by trained experts authorized by Miltenyi Biotec. An improperly serviced instrument may lead to severe personal injury or death by exposure to mechanical hazards, hazardous optical radiation, and the risk of electric shock. Unplug the power cord before cleaning the instrument. Do not allow fluids or objects to enter the interior of the instrument. Always wear protective gloves and eyewear in order to protect against potential biohazard exposure. Dispose of contaminated tissues and swabs according to their biological security level.

8.1 Decontaminating the MACSQuant® Instrument

If the instrument has been used to process biohazardous samples, the fluidics system shall be decontaminated by running the standard shutdown sequence. It is recommended to use the appropriate disinfectant for the potential pathogen. Refer to chapter 9.5 **Excessive debris is present in acquisition** if you would like to perform the **Clean** program with MACS Bleach Solution. Use tissues or swabs to decontaminate surfaces. The uptake needle of the robotic arm and the surface of the instrument shall be decontaminated upon contact with biohazardous samples.

8.2 Exchanging fluid bottles

WARNING! Do not unscrew the lids until the bottles are placed in the bottle holder to avoid spilling on the instrument and risk of electric shock. When working with biohazardous samples, it is recommended to fill the waste bottle with 100 mL of disinfectant (e.g. MACS Bleach Solution) before use. Empty bottles from the Running Buffer, Washing Solution, or Storage Solution may be used as waste bottles only with the appropriate safety label on it. For proper disposal, please follow local regulations and carefully read the chapter **Important safety information**.

- 1 Install one fluid bottle at a time. Remove the bottle to be exchanged from the bottle holder.
- 2 Place a new, sealed bottle into the orange bottle holder. Please note the corresponding color-coding (see table 2.1).
- 3 Unscrew the lid of the new fluid bottle.
- 4 Remove the bottle closure from the old bottle and install it on the new bottle. Do not disconnect the color-coded tubing and sensor cable from the bottle closure. Take care not to invert the bottle closures as this could cause liquid to enter the Hydrophobic Air Filter.

8.3 Exchanging the MACSQuant® Column

The MACSQuant® Column should be exchanged every three months or replaced with the column substitute after three months. Exchange the column as outlined in chapter 2.4.4 **Installation of the MACSQuant Column**.

8.4 Exchanging the Hydrophobic Air Filters

Hydrophobic Air Filters are attached to the bottle closures to vent the liquid bottles. To avoid clogging of the filters and to prevent contamination of the liquids, the Hydrophobic Air Filters should be exchanged if they come into contact with liquid. They also should be exchanged once a year to avoid clogging through the deposition of dust. Please use only the Hydrophobic Air Filters from Miltenyi Biotec.

8.5 Cleaning the washing station

- 1 Remove the Single Tube Rack if needed.
- 2 Open the cover of the washing station to the left side.
WARNING! Potentially contaminated liquid may spill out of the orifice of the washing station and the tubing.
- 3 Clean the washing station by wiping it with tissue and an appropriate disinfectant, e.g. 70% ethanol, isopropyl alcohol, or 10% bleach. Finally, rinse with distilled water.

8.6 Cleaning the uptake needle

The uptake needle of the robotic arm should be cleaned regularly in order to prevent contamination or blockages. Clean the uptake needle by wiping it with tissue and an appropriate disinfectant, e.g. 70% ethanol, isopropyl alcohol, or 10% bleach. Finally, rinse with distilled water. Please note that if you move the robotic needle arm up and down to get access to the entire surface of the uptake needle, the robotic needle arm has to be calibrated again. Please refer to chapter 5.1.1 **Calibration of the robotic needle arm** for more information.

8.7 Long-term storage of the MACSQuant® Instrument

The MACSQuant® Instrument should be run every two weeks in order to prevent blockages. To prevent this, every two weeks turn on the instrument and switch it into acquisition mode. This will initiate the priming of all of the fluidics. After the instrument has been primed and all system checks have been performed successfully, shut down the instrument manually as explained in chapter 3 **Switching on/off the MACSQuant Instrument**. Repeat this procedure once every two weeks.

9

Troubleshooting

Read the chapter **Important safety information**. Do not service the equipment yourself. Service shall be carried out only by trained experts authorized by Miltenyi Biotec. An improperly serviced instrument may lead to severe personal injury or death by exposure to mechanical hazards, hazardous optical radiation, and the risk of electric shock. Unplug the power cord before cleaning the instrument. Do not allow fluids or objects to enter the interior of the instrument. Always wear protective gloves and eyewear in order to protect against potential biohazard exposure. Dispose of contaminated tissues and swabs according to their biological security level.

The MACSQuant Instrument automatically analyzes the functionality of its hardware components during the instrument initialization procedure. Please contact Miltenyi Biotec Technical Support if an error message appears. This section addresses problems that are not indicated by a warning or error message that might occur during measurement, cell separation, or rinsing programs.

9.1 Column leakage

- 1 If a freshly installed MACSQuant®Column shows signs of leakage, check if the column is installed properly. Check also that the column is not broken or cracked. The column should be inserted correctly into the column connector and fastened to the point of resistance. If this is not the case, loosen the column connector, insert the column correctly, and tighten the connector.
- 2 If the column leakage persists, unscrew the column and check if the connectors of the columns are damaged. If this is the case, exchange the leaking column with a new MACSQuant Column.
- 3 If the problem still persists, contact Miltenyi Biotec Technical Support.

9.2 Syringe pump leakage

Depending on the level of use, the syringe pumps (see figure 1.1) should be cleaned every 1 to 3 months. Correct overnight and long-term storage assures that no salt deposits accumulate in the syringe pumps, which could cause wear of the pump seal. The syringe pump should not run dry at any time. This damages the pump seal. In case of leakage of the syringe pump turn off the instrument and check the connection between the syringe pump and the valve. If it is loose, tighten it. If the leakage persists, contact Miltenyi Biotec Technical Support to replace either the pump syringe (# 130-090-339) or the pump seal (# 130-022-101).

9.3 MACS® MiniSampler does not move properly

Check whether the MACS® MiniSampler guide is inserted properly into the corresponding slot on the MACSQuant Instrument and ensure that the cable connecting the MiniSampler to the back of the instrument is not obstructed in any way. See also chapter 2.4.3 **Installation of the MACS MiniSampler**.

9.4 Air bubbles during measurement or failure of event acquisition

Check for leakage within the fluidics and tubing system and perform a rinse program. Also, ensure that the sheath particle filter is devoid of air. Ensure that the needle calibration was performed correctly. Furthermore, check that the tubes connected to the bottle closures are tightly sealed. If the problem persists, contact Miltenyi Biotec Technical Support.

9.5 Excessive debris is present in acquisition

- 1 Check all fluidics to ensure that all tubes are connected properly. Ensure that the sheath particle filter is devoid of air.
- 2 Rinse or clean the system. Push the **Rinse** button in the **Instrument status bar**. Perform this twice.
- 3 If this does not work, try to clean the system by using the **Clean** and/or **Flush** programs, which can be accessed by a right-click on the **Rinse** button. Press the **Clean** button to perform a system clean with MACS Bleach Solution. Follow the instructions in the dialog box that appears and place the required bleach solution in the Single Tube Rack. The Single Tube Rack must be attached and selected in the **Experiment** tab in order to perform a **Clean** program. The **Flush** program will perform a system backflush. This process will take approximately 16 minutes.

9.6 Touchscreen remains dark

If the touchscreen remains dark, check if the power cord is plugged in and the power supply is switched on. If the instrument is switched on and the touchscreen remains dark, switch off the instrument, wait five seconds, and switch on again. If the MACSQuant Instrument still does not initialize, contact Miltenyi Biotec Technical Support.

10

Technical data and specifications

10.1 Technical data and specifications of the MACSQuant[®] Instrument

The MACSQuant[®] Instrument is labeled as a protection class I device and must be plugged into a grounded power outlet.

WARNING! Read the chapter **Important safety information** before installation and use. The MACSQuant Instrument is not specified for use in a cold room.

Please note that for USA and Canada the main power supply cord and plug of the instrument shall comply with following specifications: UL listed and KAM cord, minimum type SJ, minimum 18 AWG, 3 conductors. Rated for a minimum temperature of 60 °C. Provided with grounding-type (NEMA 5-15P) attachment plug, rated 125 VAC, 10 A. Opposite end terminates in IEC 320 style connector, rated 125 VAC, 10 A.

Parameter	Specification
Models	MACSQuant Analyzer (# 130-092-197) MACSQuant Analyzer 10 (# 130-096-343) MACSQuant VYB (# 130-096-116)
Color	Green/orange
Footprint	605 mm × 400 mm (w × d)
Footprint with MACS MiniSampler	605 mm × 500 mm (w × d)
Height	392.5–520 mm (adjustable touchscreen)
Weight	50 kg
Input voltage	100–240 VAC, 50/60 Hz
Power consumption	350 VA
Fuses	2 × 8 AT, 250 V
RS232 Interface (labeled “RS232/AUX”) Not in use	Pin 1, 4, 6–9: NC Pin 2: RXD Pin 3: TXD Pin 5: GND
AC Output (labeled “Bottle Sensor”)	Pin 1–5: Input Pin 6–15: GND Pin 9–13: Input/Output 5 VAC/10 kΩ
RS232 Interface + DC-Output (labeled “RS232/BCR”)	Pin 1: Input/Output Pin 2: TXD Pin 3: RXD Pin 4, 6: NC Pin 5: GND Pin 7, 8: Shorted Pin 9: 5 VDC/2.6 A
CAN Bus + DC-Output (labeled “External CAN”)	Pin 1, 4, 8: NC Pin 2: CAN-L Pin 3, 6: GND Pin 5, 9: 24 VDC/1.85 A Pin 7: CAN-H
RS232 Interface (labeled “COM 1”)	Pin 1: DCD Pin 2: RxD Pin 3: TxD Pin 4: DTR Pin 5: GND Pin 6: DSR Pin 7: RTS Pin 8: CTS Pin 9: RI
4 × USB 2.0 ports (labeled “USB 1-4”)	Pin 1: +5 V Pin 2: - Pin 3: + Pin 4: GND

Table 10.1: Technical specifications of the MACSQuant Instrument.

Parameter	Specification
VGA port (labeled "VGA")	Video graphics array analog interface supports up to 2048 × 1536 @75 Hz Pin 1: RED Pin 2: GREEN Pin 3: BLUE Pin 4: ID2/RES Pin 5-8, 10: GND Pin 9: KEY/PWR Pin 11: ID0/RES Pin 12: ID1/SDA Pin 13: HSync Pin 14: VSync Pin 15: ID3/SCL
2 × HDMI ports	2 × High-Definition Multimedia Interfaces support up to 1920 × 1200 @60 Hz Pins 1-3: TMDS Data2 Pins 4-6: TMDS Data1 Pins 7-9: TMDS Data0 Pins 10-12: TMDS Clock Pin 13: CEC Pin 14: Reserved Pin 15: SCL Pin 16: SDA Pin 17: GND Pin 18: +5 V/50 mA Pin 19: Hot Plug detect
2 × Ethernet ports (GbE) (labeled "LAN 1" and "LAN 2")	Pin 1: TXD+ Pin 2: TXD- Pin 3: RXD+ Pins 4, 5, 7, 8: NC Pin 6: RXD-
Audio ports	Green: headphones Pink: microphone
CF Card Slot	The inserted card contains the operating system of the MACSQuant Instrument (Microsoft Windows XP with MACSQuantify Software, read-only mode)
RAM	2 GB DDR3 (SO-DIMM, 1066 MHz)
HDD	2.5" automotive & extended temperature range
Monitor	12" TFT touchscreen
Working temperature	+20 °C to +25 °C
Humidity	0% to 85% relative humidity, non-condensing
Altitude	Max. 2000 m
Emission sound pressure level at workstation	<70 dB(A)

Table 10.2: Technical specifications of the MACSQuant Instrument continued.

Parameter	Specification
Minimal sample volume	25 µL (recommended)
Minimum sample particle size	1 µm
Sample flow rate	25–100 µL/min plus automated flow rate to maintain 500, 1000, or 2000 events/second
Sheath fluid consumption	2–10 mL/min
Maximal event rate	Up to 10 ⁴ events/second
Lasers MACSQuant Analyzers	405 nm, 488 nm, 635 nm
Lasers MACSQuant VYB	405 nm, 488 nm, 561 nm
Fluorescence sensitivity	MESFs (CV < 5%); FITC < 200; PE < 100; APC < 150
Absolute cell count performance	Accuracy +/- 5%; reproducibility (CV) <5%
Sample carryover	Typical value of 0.01%

Table 10.3: Operation details of the MACSQuant Instrument.

Model	Channel	Laser [nm]	Filter
MACSQuant Analyzer	FSC	488	488/10
	SSC	488	488/10
	V1	405	450/50
	B1	488	525/50
	B2	488	585/40
	B3	488	655-730
	B4	488	750LP
	R1	638	655-730
	R2	638	750LP
MACSQuant Analyzer 10	FSC	488	488/10
	SSC	488	488/10
	V1	405	450/50
	V2	405	525/50
	B1	488	525/50
	B2	488	585/40
	B3	488	655-730
	B4	488	750LP
	R1	638	655-730
R2	638	750LP	
MACSQuant VYB	FSC	561	561/10
	SSC	561	561/10
	V1	405	450/50
	V2	405	525/50
	B1	488	525/50
	B2	488	614/50
	Y1	561	586/15
	Y2	561	615/20
	Y3	561	661/20
Y4	561	750LP	

Table 10.4: Overview of the different MACSQuant Instrument configurations.

Supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage. Transient over voltages (voltage in excess of the normal operating voltage) present on the mains supply: Category II. The instrument is suitable for rated pollution degree 2.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications of the equipment unless expressly approved by Miltenyi Biotec may void your authority to operate the equipment pursuant to FCC 47 CFR.

The MACSQuant Instrument has been investigated by Underwriter Laboratories in accordance with the standards

- UL 61010-1
- CAN/CSA-C22.2 No. 61010-1
- IEC 61010-2-081

and complies with all essential requirements of the following European directives

- 2004/108/EC (electromagnetic compatibility)
- 2006/42/EC (machinery)

and is in conformity with the following harmonized standards

- EN 61010-1
- EN 61010-2-081
- EN 61326-1
- EN 62471.

For other safety considerations, refer to the product label, or visit **www.miltenyibiotec.com**.

Design and specifications are subject to change without notice.

10.2 Technical data and specifications of the MACS® MiniSampler

The MACS® MiniSampler is labeled as a protection class III device and must only be plugged into the connector labeled with External CAN of the MACSQuant Instrument.

WARNING! Read the chapter **Important safety information** before installation and use. The MACS MiniSampler is not specified for use in the cold room.

Parameter	Specification
Size without lid	182 mm × 148 mm × 47 mm
Size with lid	280 mm × 153 mm × 172 mm
Weight	1.5 kg
Input voltage	24 VDC
Current	0.8 A
Sub D9 interface with shielding	Pins 1, 4, and 8: NC Pin 2: CAN-L Pin 3, 6: GND Pin 5, 9: +24 VDC/2A Pin 7: CAN-H
Working temperature	+20 °C to +25 °C
Humidity	0% to 85% relative humidity, non-condensing
Altitude	Max. 2000 m

Table 10.5: Technical data and pin assignment of the MACS MiniSampler.

The MACS MiniSampler has been investigated by Underwriter Laboratories in accordance with the standards

- UL 61010-1
- CAN/CSA-C22.2 No. 61010-1

and complies with all essential requirements of the following European directives

- 2004/108/EC (electromagnetic compatibility)
 - 2006/42/EC (machinery)
- and is in conformity with the following harmonized standards
- EN 61326-1.

11

Technical Support

For technical support, please contact your local Miltenyi Biotec representative or the Technical Support team at Miltenyi Biotec headquarters:

Miltenyi Biotec GmbH

Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-830
Fax +49 2204 8306-89
macstec@miltenyibiotec.de

Visit **www.miltenyibiotec.com/local**
to find your nearest Miltenyi Biotec contact.

12

Limited warranty

Except as stated in a specific warranty statement, which may accompany this product, or unless otherwise agreed in writing by an authorized representative of Miltenyi Biotec, Miltenyi Biotec's warranty, if any, with respect to this product is subject to the terms and conditions of sale of the company within the Miltenyi Biotec group which supplied the product. Terms and conditions of sale may vary by country and region. Nothing in this document should be construed as constituting an additional warranty.

Miltenyi Biotec's warranty for this product only covers product issues caused by defects in material or workmanship during ordinary use; it does not cover product issues caused by any other reason, including but not limited to product issues due to use in a manner other than specifically described in this manual, for example: inappropriate or improper use; incorrect assembly or installation by an operator or a third party; reasonable wear-and-tear; negligent or incorrect handling, servicing, or maintenance; non-adherence to the operating instructions; unauthorized modification of or to any part of this product; or use of inappropriate accessories or work materials.

Miltenyi Biotec's warranty does not cover products sold AS IS or WITH ALL FAULTS, or which had its serial number defaced, altered or removed, or any consumables, or parts identified as being supplied by a third party. Miltenyi Biotec must be informed immediately, if a claim is made under such warranty. If a material or manufacturing defect occurs within the warranty period, Miltenyi Biotec will take the appropriate steps to restore the full usability of your product.

Limitation on damages: Miltenyi Biotec shall not be liable for any incidental or consequential damages for breach of any express or implied warranty or condition on this product.

Some states or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations or exclusions may not apply to you. This warranty statement gives you specific legal rights and you may have other rights, which may vary from country to country or jurisdiction to jurisdiction.

13

Glossary

Air filter connector	Luer-to-thread connector for attaching the Hydrophobic Air Filters to the threaded bottle closure vent.
Bottle closure	Vented screw-on closure with fluid uptake tubes. The bottle closures contain fluid sensors and are equipped with sensor cable connectors. Each bottle closure is color coded
Cell enrichment	Process of cell enrichment is accomplished by the MACS Enrichment Unit and the MACSQuant Column. Target cells labeled with MACS MicroBeads are enriched using MACS Technology.
Chill Racks	Different acrylic tube racks are available with the instrument and are designed for optimal positioning of sample tubes. They contain a coolant allowing racks to be pre-cooled in the refrigerator for subsequent cooling of the cells during analysis.
Column connector	Luer-to-thread connector holding the MACSQuant Column.
Column substitute	Column devoid of any paramagnetic particles. This column substitute can be used when enrichments are not being performed or for long-term storage and shipment. The column substitute cannot be used for cell enrichment prior to analysis.

Fluid bottle	1.5 L bottle that holds the system fluids for operational use of the MACSQuant Instrument. Fluid sensors monitor the fluid levels in the bottles for Running Buffer, Washing Solution, Storage Solution, and waste via electrolyte conductivity.
Fluid sensor cable	Cable connecting the fluid sensor to the MACSQuant Instrument. The fluid sensor cables are color-coded: blue for Running Buffer, green for Washing Solution, black for Storage Solution, and red for waste.
Fluid sensors	The fluid sensors measure electrolyte conductivity and is integrated into the closures of fluid bottles and the waste bottle.
Front panel	The front panel opens sideways, giving access to the MACSQuant Instrument Column, pumps, valves, washing station, and tubings.
Hydrophobic Air Filter	Hydrophobic Air Filter (0.2 μm) attached to the bottle closure. Used to vent the bottle and simultaneously prevent contaminants from entering or escaping from the bottle.
MACS Enrichment Unit	Black cover surrounds the area enclosing the magnet. The magnet cover is located in the center of the fluidics system and has a slot for the insertion of the MACSQuant Column.
MACS MicroBeads	Superparamagnetic particles conjugated to antibodies for magnetic labeling of cells or biomolecules.
MACS Technology	Technology for immunomagnetic labeling and subsequent separation of cells or biomolecules in a high-gradient magnetic field.
MACSQuant Column	Separation column specifically designed for the MACSQuant Instrument.
MACSQuant Instrument	An automated flow cytometer that is also referred to as device or instrument.
MACSQuant Running Buffer	Sterile and ready-to-use buffer for flow cytometry, cell enrichment, and washing programs. The tubing connector is color-coded blue.
MACSQuant Storage Solution	Sterile and ready-to-use solution for overnight and longterm storage of the MACSQuant Instrument. The tubing connector is color-coded black.

MACSQuant Washing Solution	Sterile and ready-to-use solution for washing and special cleaning programs. The tubing connector is color-coded green.
Negative fraction	Sample fraction containing the cells not magnetically labeled. During MACS Separation, these cells are not retained on the column and pass directly through the column while the column is in the presence of the MACS Magnet.
Positive fraction	Sample fraction containing the cells labeled with MACS MicroBeads after MACS Separation. These cells are retained on the column while the column is placed in the magnetic field. The cells are eluted from the column after the column has been removed from the magnet.
Sheath particle filter	Sheath particle filter located on the back of the MACSQuant Instrument.
Syringe pumps	Computer-controlled, high precision syringe pumps with Teflon® seal plunger that drive fluids through the MACSQuant Instrument fluidics system.
Touchscreen	High resolution TFT color touchscreen located on top of the MACSQuant Instrument. The touchscreen is used to operate and monitor the instrument via on-screen menus.
Tubing connector	Threaded plastic connector with a square nut used to connect the tubings to the bottle closures, the column, the pump, or valves.
Tubing system	Permanent set of Teflon® tubing through which fluid circulates in the MACSQuant Instrument fluidics system
Uptake needle	The robotic needle arm carries a needle designed for computer-controlled sample uptake. This uptake needle is calibrated and self-cleaned. It is designed to move in the y- and z-direction.
Waste bottle	Bottle for waste fluid. The closure is equipped with a fluid sensor. The closure, the fluid sensor cable, and the tubing connector are color-coded red.

Miltenyi Biotec

Germany/Austria/ Switzerland

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

USA/Canada

Miltenyi Biotec Inc.
2303 Lindbergh Street
Auburn, CA 95602, USA
Phone 800 FOR MACS
Phone +1 530 888 8871
Fax +1 530 888 8925
macs@miltenyibiotec.com

www.miltenyibiotec.com

Australia

Miltenyi Biotec
Australia Pty. Ltd.
Unit 16A, 2 Eden Park Drive
North Ryde, NSW 2113,
Australia
Phone +61 2 8877 7400
Fax +61 2 9889 5044
macs@miltenyibiotec.com.au

Benelux

Miltenyi Biotec B.V.
Schipholweg 68 H
2316 XE Leiden
The Netherlands
macs@miltenyibiotec.nl

Customer service The Netherlands

Phone 0800 4020120
Fax 0800 4020100

Customer service Belgium

Phone 0800 94016
Fax 0800 99626

Customer service Luxembourg

Phone 800 24971
Fax 800 24984

China

Miltenyi Biotec GmbH
Shanghai Office
Rm. 2309
No. 319 Xianxia Rd.
Shanghai 200051, P.R. China
Phone +86 21 62351005
Fax +86 21 62350953
macs@miltenyibiotec.com.cn

France

Miltenyi Biotec SAS
10 rue Mercoeur
75011 Paris, France
Phone +33 1 56 98 16 16
Fax +33 1 56 98 16 17
macs@miltenyibiotec.fr

Italy

Miltenyi Biotec S.r.l.
Via Persicetana, 2/D
40012 Calderara di Reno (BO)
Italy
Phone +39 051 6 460 411
Fax +39 051 6 460 499
macs@miltenyibiotec.it

Japan

Miltenyi Biotec K.K.
Nittsu-Eitai Building 5F
16-10 Fuyuki, Koto-ku,
Tokyo 135-0041, Japan
Phone +81 3 5646 8910
Fax +81 3 5646 8911
macs@miltenyibiotec.jp

Nordics and Baltics

Miltenyi Biotec Norden AB
Scheelevägen 17
223 70 Lund
Sweden
macs@miltenyibiotec.se

Customer service Sweden

Phone 0200-111 800
Fax 046-280 72 99

Customer service Denmark

Phone 80 20 30 10
Fax +46 46 280 72 99

Customer service Norway, Finland, Iceland, and Baltic countries

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Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS, MACSQuant, and MACSQuantify are registered trademarks or trademarks of Miltenyi Biotec GmbH. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only. Copyright © 2013 Miltenyi Biotec GmbH. All rights reserved.

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Asia Pacific Pte Ltd.
Miltenyi Biotec Ltd.