

NT-proBNP Validations

- Validation includes:
- % detectability data for Supplier B's Human CVD3 multiplex panel

Suppliers

- A: Bio-Rad
- B: Linco (done)
- C: Upstate
- D: Biosource
- E: RnD Systems

NT-proBNP - % Detectability

- Using supplier B's Human CVD3 panel, % detectability of NT-proBNP in serum and plasma from healthy donors:
 - EDTA plasma (n=10) 60%
 - Citrated plasma (n=10) 0%
 - Serum (n=10) 30%
- (Note that for each sample type a different set of donors were used.)