

VEGF Validation

- Validation includes:
- % detectability data for Supplier B's Human CVD3 multiplex panels

Suppliers

- A: Bio-Rad
- B: Linco (done)
- C: Upstate
- D: Biosource
- E: RnD Systems

VEGF - % Detectability

- Using supplier B's Human CVD3 panel, % detectability of VEGF in serum and plasma from healthy donors:
 - EDTA plasma (n=10) 40%
 - Citrated plasma (n=10) 0%
 - Serum (n=10) 10%
- (Note that for each sample type a different set of donors were used.)